



Step-by-Step Clinical Guide

Appointment 1 - Initial Scan

Full intraoral scan (edentulous arches + bite).

Clean, dry, detail-focused scan.

Appointment 2 - Try-in

3D-printed denture try-in for occlusion and esthetic adjustments.

Requires secondary Rx.

Appointment 3 - Final Seat

Delivery of finalized milled or printed digital denture.

Make minor adjustments chairside.

Digital Denture Scan Requirements
Preparation: Clean and Dry Field: Use air to fully dry tissues; saliva will distort edentulous scans.

Scanning Steps (Edentulous/No Teeth)

Table with 3 columns: Area, What to Capture, How. Rows include Start Point, Upper Arch, Lower Arch, and Bite Registration.